



# Nano-Spray-Ionization Source




- The NSI source can be adjusted to fit all Thermo instruments (LTQ, Orbitrap, TSQ, Fusion)
- It offers all the benefits and features of commercial NSI sources:
  - Lower flow rates: no drying gas or thermal heating is required and increased sensitivity compared to higher flow ESI
  - Higher tolerance to a wide variety of liquid compositions than conventional ESI
- Additionally this NSI source offers significant improvements over most commercial NSI sources:
  - Great flexibility: it will not only accommodate in house packed columns of all length but also commercial columns, as well as configurations with traps, spray tips etc.
  - The source can be customized to your specific needs: e.g. mounting of cameras, LED lights or anything else you can think of.....
  - Additionally this NSI source can be used in conjunction with the Thermo NSI probes (static , dynamic and packed tip probe)
  - The open design allows for leaks to be detected immediately
  - Significantly lower cost
  - The NSI platform can also be used with the adaptor ring of the Thermo NSI source, eliminating the need to machine the adaptor bracket
  - All the parts are commercially available , except for the adaptor system, which is custom machined, e.g. by emachineshop.com
- The UWPR is currently using the NSI source on our Thermo instruments
- The plans for this NSI source are available for free, use at your own risk.... 😊
  
- This document and the xls file with all the part numbers should help you build your own source
- Note this is just a basic version, you can modify this source to fit your needs


This is what we are going to build:




# UWPR NSI Source



# UWPR NSI Source



# UWPR NSI Source




# NSI Source: step 1



**Note: you need two sets per instrument (one for each pin)**

1. Attach base to post holder using a ¼-20 screw
2. Slide the bushings into the post holder in the order shown above
3. Replace the thumb screw with a ¼-20 set screw (optional)
4. Lightly tighten the set screw (or thumb screw) to hold flanged bushing
5. Slide the assembled post holders over the pins on the instrument

## NSI Source: step 2



5. Slide the two post holders over the two pins on the mass spectrometer
6. While attached to the mass spec use two 1/4-20 screws to attach the construction rail to the two bases  
Make sure the construction rail is level and tighten the screws

CR-1

Construction Rail, 1 x 5 in.

(Newport)



# NSI Source: step 3

Use one #8-32 x 1/4" with a #8-32 washer and nut and one #8-32 x 3/8" with a nut, to attach 2.5" rail to rail carrier




Newport 2.5" rail (9022NF)

Newport rail carrier (9025NF)

Replace #8-32 setscrew with Thumbscrew (91035A430)

Use two button head screws you received with the NewFocus rail, to attach rail to the breadboard

Use two 1/4"-20 x 3/4" screws to attach breadboard to CR-1



Newport 12" rail (9020NF) to mount  
Newport rail carrier (9025NF) with attached 2.5" rail (9022NF)


- |           |  |                |
|-----------|--|----------------|
| 9022NF    | e-z-Trac™ Rail, 2.5 in. Length                         | (Newport)      |
| 9025NF    | e-z-Trac™ Post Carrier, 0.75 in. Width                 | (Newport)      |
| 91035A430 | 18-8 Stainless Steel Knurled Head Captive Panel Screw  | (McMasterCarr) |
| 9020NF    | e-z-Trac™ Rail, 12.0 in. Length                        | (Newport)      |
| MB810     | Aluminum Breadboard, 8" x 10" x 1/2", 1/4"-20 Threaded | (Thorlabs)     |

# NSI Source: step 4




Mount Newport XYZ stage, part# 9065-XYZ  
To 2.5" rail using set screws provided with stage

Attach angle plate, part# 9101 to XYZ stage using  
two 1/4"-20 x 3/8" screws



set screws




UWPR NSI source

9065-XYZ  
AJS100-0.5K-NL  
9101NF

Pint-Sized XYZ-Translation Stage  
High Precision Large Knob Adjustment Screw, 12.7 mm Travel, 100 TPI, No Lock  
Angle Plate, 1/4-20 Thru Slots, 8-32 and 1/4-20 Threads

(Newport)  
(Newport)  
(Newport)

# NSI Source: step 5



RLA1200  
RC1  
VC1  
T711-250  
94323A598

Imperial Dovetail Optical Rail, 12"  
Rail Carrier, Counterbored Hole 1" x 1"  
Small V-Clamp with PM3 Clamping Arm  
Post Mountable Fiber Clamp, 250mm  
Nylon Raised-Head Thumb Screws  
8-32 Thread Size, 1" Long

(Thorlabs)  
(Thorlabs)  
(Thorlabs)  
(Thorlabs)  
(McMasterCarr)

- ▶ Adjustable Force Magnetic Clamp for 250µm Jacketed Fiber
- ▶ Clamping Arm Swings Clear for Easy Loading
- ▶ Rubber Pad Provides Excellent Holding Force

Magnetically Coupled Swing Arm Makes Loading/Unloading Fiber Simple

T711-250

Soft Rubber Pad Securely Clamps 0250µm Fiber In Machined V-Groove

T711-250

Related Products

#8-32 (M4) Tap

Multimode Fiber

Bare Fiber Terminator

Fiber Launch Systems



# NSI Source: step 6a

## To bypass interlock on ThermoElectron LTQ, TSQ, QE when using alternate nanospray source (NSI)

**Note:** Jumper and resistor can either be attached directly to the connector on the front of the LTQ or via a DB15 male connector (easier to change between different sources).


Jumper to bypass interlock between socket 9 and 10

10 K ohm resistor (2% var.) for NSI source recognition between socket 7 and 8



### OR use DB15 male connector:

solder the jumper and resistor to the corresponding pins



DB15 male connector backside view



Solder Cup D-Sub Connector, DB15 Male



510-2002  
510-2445  
70201054

Emerson Network Power, D-Sub Plug; Thermoplastic; Plug; Copper Alloy; Solder Termination  
Emerson Network Power, Hood; D-Sub; Chrome; Metalized Plastic; Steel; UL 94 VO; RoHS Compliant  
Resistor; Metal Film; Res 10 Kilohms; Pwr-Rtg 0.25 W; Tol 2%; Axial; Epoxy

(Allied Electronics)  
(Allied Electronics)  
(Allied Electronics)

# NSI Source: step 6b

To bypass interlock on ThermoElectron Fusion or Quantiva when using alternate nanospray source (NSI)



Male Plug HARTING 09200102612: Connect 10 kΩ resistor between pin 6 and 10 and jumper wire between pin 4 and 9  
Optionally a ground wire can be attached to pin 7

In case you wonder why I show it upside down... that's because looking at the instrument from the front, that is how it will plug into the instrument


70104585  
70201054

Insert; 250 V; 16 A; Male; 10; RoHS Compliant; Han A Series; 10 A Product Size Mfr. Part#: 09200102612  
Resistor; Metal Film; Res 10 Kilohms; Pwr-Rtg 0.25 W; Tol 2%; Axial; Epoxy

(Allied Electronics)  
(Allied Electronics)

# NSI Source: step 7


**IMPORTANT: cover exposed metal with heat shrink polyolefin tubing**



10334-02	SCID series, threaded plastic flange receptacle assembly	(Connectronics Corp)
Or		
00004-95039	Connector, Feed thru, 15KVdc, 7.5A, Connectronics	(Unity Lab Services)
066-031-14	High voltage wire, #22 AWG, .180 O.D. silicone /rubber insul.	(Connectronics Corp)
512-2474	20DF PIN CONT. LP	(Allied Electronics)
512-2478	20DF SOCKET CONT. L.P.	(Allied Electronics)
AA43288BU	Platinum wire; 99.95% ; Alfa Aesar; 0.5mm dia; 25cm; 4.21g/m	(Fisher)
P-777	Microcross, .025" OD tubing sleeves, .006" THRU HOLE, PEEK	(Upchurch)
P-875	Microtee , w/ mount. whole, for .025" OD tubing sleeves, PEEK™	(Upchurch)
F-172x	Ferrule, for .025" sleeves/P-416 NUT, PEEK™, BLACK (10 PK)	(Upchurch)
F-185x	Tubing sleeve, 395µm (.015") ID x .025", PEEK™, GREEN (10 PK)	(Upchurch)
P-116	Plug, for micro fittings, PEEK™, BLACK	(Upchurch)

# NSI Source: step 8


Note on the Fusion and Quantiva there is a inter-lock that needs to be pushed in  
I used the miniature posts and a 0.5" 4-40 set screw



Adjust set screw to make sure the interlock is pressed in



MS3R



Set screw 4-40 0.5"



MSA25



ER90B


MS3R  
MSA25  
ER90B  
Set screw 4-40 0.5"

Mini Series Mounting Posts, 6mm Diameter, 3" Long  
Thread Adapter, 1/4"-20 to #4-40  
Mini-Post Right Angle Adapter

(Thorlabs)  
(Thorlabs)  
(Thorlabs)

# NSI Source: step 9

Note on the QE there is a inter-lock that needs to be pushed in  
I used the miniature posts and a 1" 4-40 set screw



Adjust set screw to make sure the interlock is pressed in



MS3R  
MS2R  
MSA25  
ER90B  
4-40 screw 1in long

Mini Series Mounting Posts, 6mm Diameter, 3" Long  
Mini Series Mounting Posts, 6mm Diameter, 2" Long  
Thread Adapter, 1/4"-20 to #4-40  
Mini-Post Right Angle Adapter

(Thorlabs)  
(Thorlabs)  
(Thorlabs)  
(Thorlabs )